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1. Drosten C et al. Clinical features and virological analysis of a case of Middle East respiratory syndrome coronavirus infection. *Lancet* 2013; published online June 17. [http://dx.doi.org/10.1016/S1473-3099\(13\)70154-3](http://dx.doi.org/10.1016/S1473-3099(13)70154-3)
2. Guery B et al. Clinical features and viral diagnosis of two cases of infection with Middle

Eastrespiratory syndrome coronavirus: a report of nosocomial transmission. *Lancet* 2013; published online May 29. [http://dx.doi.org/10.1016/S0140-6736\(13\)60982-4](http://dx.doi.org/10.1016/S0140-6736(13)60982-4)

3. Memish ZA et al. Respiratory tract samples, viral load, and genome fraction yield in patients with Middle East respiratory syndrome. *J Infect Dis* 2014; published online May 15. [http://J Infect Dis. \(2014\)doi:10.1093/infdis/jiu292](http://journals.pubs.lit.edu/doi/10.1093/infdis/jiu292)

4. Kapoor M et al. Clinical and laboratory findings of the first imported case of Middle East respiratory syndrome coronavirus (MERS - CoV) into the United States. *Clin Infect Dis* 2014; published online August 6. <http://cid.oxfordjournals.org/content/early/2014/08/05/cid.ciu635.long>

5. Assiri A et al. Epidemiological, demographic, and clinical characteristics of 47 cases of Middle East respiratory syndrome coronavirus disease from Saudi Arabia: a descriptive study. *Lancet* 2013; published online July 26. [http://dx.doi.org/10.1016/S1473-3099\(13\)70204-4](http://dx.doi.org/10.1016/S1473-3099(13)70204-4)

6. Kraaij – Dirkzwager M et al. Middle East respiratory syndrome coronavirus (MERS - CoV) infections in two returning travellers in the Netherlands, May 2014. *Eurosurveillance* 2014; 19:20817 <http://www.eurosurveillance.org/ViewArticle.aspx?ArticleId=20817>

7. Abroug F et al. Family Cluster of Middle East Respiratory Syndrome Coronavirus Infections, Tunisia, 2013. *Emerg Infect Dis* 2014 Sep. <http://dx.doi.org/10.3201/eid2009.140378>

8. Corman VM et al. Detection of a novel human coronavirus by real-time reverse-transcriptase polymerase chain reaction. *Eurosurveillance* 2012; 17:20285. <http://www.eurosurveillance.org/ViewArticle.aspx?ArticleId=20285>

9. Corman VM et al. Assays for laboratory confirmation of novel human coronavirus (hCoV-EMC) infections. *Eurosurveillance* 2012; 17:20334. <http://www.eurosurveillance.org/ViewArticle.aspx?ArticleId=20334>

10. Novel Coronavirus 2012 Real-Time RT-PCR Assay, Centers for Disease Control and Prevention. Available at: <http://www.fda.gov/downloads/MedicalDevices/Safety/EmergencySituations/UCM355572.pdf> Accessed 25 July 2014.

11. Buchholz U et al. Contact investigation of a case of human coronavirus infection treated in a German hospital, October-November 2012. *Eurosurveillance* 2013 ;18:20406. <http://www.eurosurveillance.org/ViewArticle.aspx?ArticleId=20406>

12. Reusken C et al. Specific serology for emerging human coronaviruses by protein microarray. *Eurosurveillance* 2013;18:20441. <http://www.eurosurveillance.org/ViewArticle.aspx?ArticleId=20441>

13. Perera RA, et al. Seroepidemiology for MERS coronavirus using microneutralisation and pseudoparticle virus neutralisation assays reveal a high prevalence of antibody in dromedary camels

in Egypt, June 2013. *Eurosurveillance* 2013;18:20574.
<http://www.eurosurveillance.org/ViewArticle.aspx?ArticleId=20574>

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